Link do produktu: <http://www.novazym.sklep.pl/novabeads-total-rna-purification-kit-100-preparations-3-zone-magnetic-particles-p-437.html>

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
Total RNA purification KIT magnetic particles & 3-Zone (Novazym Polska)

TISSUE PREPARATION:


1) Isolate of nucleoprotein complex by using a power homogenizer (BIOPREP-24 or equivalent) in **3-ZONE** (1 mL reagent per each 50-100 mg of tissue). The tissue sample volume should not exceed 10% of the volume of 3-ZONE reagent used for homogenization.


2) Incubate the homogenized sample for 5 minutes at 15 to 30°C to permit the complete dissociation of nucleoprotein complexes.


2) Add 0.2 mL of chloroform per 1 mL of 3-ZONE Reagent.


3) Cap sample tubes securely, and shake them vigorously by hand for 15 seconds.

4) Incubate tubes at 15 to 30°C for 2 to 3 minutes.

5) Centrifuge the samples at 12,000 × g, 15 minutes in temperature 2 to 8°C.

6) Following centrifugation, the mixture separates into a lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. The volume of the aqueous phase is about 60% of the volume of 3-ZONE Reagent used for homogenization.

RNA PURIFICATION:

1) Transfer the aqueous phase to a fresh tube, add the same volume of 70% ethanol and mixed well. Add 25 µl of magnetic beads and incubate 10 min. at room temperature.


2) Put the tubes on the magnetic rack and then remove the supernatant.

3) Transfer the tubes to an ordinary rack (without a magnet) and wash with 1 mL of 70% EtOH.

4) Repeat p.3

5) Remove residual ethanol (very carefully). Leave the open tubes for 10-15 minutes in a magnetic rack to evaporate residual ethanol.

6) Remove the tubes from the magnetic rack. Resuspend beads in 50 µl of deionised water (free of RNA). Incubate 5 minutes at 75 °C, 800 rpm **. Then, open the tubes and incubate for another 5 minutes.

7) Transfer the tubes to the magnetic rack. Supernatant, containing RNA, transfer to a new, clean and RNases free tube.

8) Measure of the RNA concentration.

* The inhibition during amplification step occurred repeat the p.3 in RNA purification protocol.

** optionally: When the elution is carried out by mean of dry block without of shaking, the eluates should be

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Fast RNA purification KIT magnetic particles & 3-Zone (Novazym Polska)

FTA card protocol

FTA SAMPLE PREPARATION:


- 1) Place a piece of FTA card (30–40 mg) in 0.5 mL reagent **3-ZONE**


- 2) Incubate sample for 5 minutes at 15 to 30°C to permit the complete dissociation of nucleoprotein complexes.

- 2) Add 0.2 mL of chloroform per 1 mL of 3-ZONE Reagent.


- 3) Cap sample tubes securely, and shake them vigorously by hand for 15 seconds.

- 4) Incubate tubes at 15 to 30°C for 2 to 3 minutes.

- 5) Centrifuge the samples at 12,000 × g 15 minutes in 2 to 8°C.

- 6) Following centrifugation, the mixture separates into a lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. The volume of the aqueous phase is about 60% of the volume of 3-ZONE Reagent used for homogenization.

RNA PURIFICATION:

- 1) Transfer the aqueous phase to a fresh tube, add 0.5 mL of 70% ethanol and mixed well. Add 25 µl of magnetic beads and incubate 10 min. at room temperature.

- 2) Put the tubes on the magnetic rack and then remove the supernatant.

- 3) Transfer the tubes to an ordinary rack (without a magnet) and wash with 1 mL of 70% EtOH.


- 4) Repeat p.3

- 5) Remove residual ethanol (very carefully). Leave the open tubes for 10–15 minutes in a magnetic rack to evaporate residual ethanol.

- 6) Remove the tubes from the magnetic rack. Resuspend beads in 50 µl of deionised water (free of RNA). Incubate 5 minutes at 75 °C, 800 rpm. Then, open the tubes and incubate for another 5 minutes.

- 7) Transfer the tubes to the magnetic rack. Supernatant, containing RNA, transfer to a new, clean and RNAses free tube.

- 8) Measure of the RNA concentration.

Novazym Polska s.c., ul. Rubież 46H 61-612 Poznań, tel. +48 0 61 610 39 10, fax +48 0 61 610 39 11, NIP 9721099648,
email info@novazym.com, www.novazym.com