

Novazym http://www.novazym.sklep.pl info@novazym.pl

Dane aktualne na dzień: 02-05-2025 08:19

Link do produktu: http://www.novazym.sklep.pl/vetpcr-nnv-96rxn-p-691.html



VetPCR? NNV (96RXN)

Numer katalogowy

VET-P003-96R

Opis produktu

Nervous necrosis virus RT-PCR Detection Kit Viral nervous necrosis is a disease found world-wide in many species of marine fish, causing an

extremely high mortality of affected larvae and juveniles. The agent of VNN has been identified as a new member of the Nodaviridae. Fish nodaviruses have been classified into four genotypes: tiger puffer nervous necrosis virus, striped jack nervous necrosis virus, barfin flounder nervous necrosis virus and red-spotted grouper nervous necrosis virus. Nerve tissue is the target tissue for nervous necrosis virus (NNV) infection, but NNV can also infect other tissues. The distribution of NNV in other tissues varies with fish species and age.

VetPCR? NNV Detection Kit is the direct detection of Nervous necrosis virus on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate, reliable technique.

Characteristics

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20?C. Time-saving and cost-effective.

Contents

Kit	Quantity 48	Quantity 96	Package
VetPCR [™] NNV RT-PCR Premixture	1	1	Vial
VetPCR™ NNV PCR Premixture	1	1	Vial
Brig [™] RT-PCR solution	1	1	Vial
Biotech™ Transcriptase solution	1	1	Vial
DNase/Rnase free water	1	1	Vial
NNV PCR Positive control	1	1	Vial
PCR Negative Control	1	1	Vial
PCR Internal Control	1	1	Vial
Mineral Oil solution	1	1	Vial
Brig™ Molecular Weight marker	1	1	Vial
RNA purification kit	50	100	Test

Interpretation of the Test Result

Expected PCR product size : 238bp



Novazym http://www.novazym.sklep.pl info@novazym.pl



Fig. 1 Result:

Lane M: Brig? Molecular Weight Marker Lane 1~2: NNV Positive samples Lane I.C.: Internal control Lane P: Positive control Lane N: Negative control