Link do produktu: <http://www.novazym.sklep.pl/ampli-lamp-reticuloendotheliosis-rev-reagent-set-96-rxn-p-487.html>

Ampli-LAMP Reticuloendotheliosis (REV) Reagent Set 96 rxn

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


1. Isothermal Master Mix (reaction buffer, $MgCl_2$, dNTPs, pyrophosphatase, GspSSD polymerase, fluorescent dye ? 6-FAM derivative).
2. Primers (mix, F3, B3, FIP, BIP, LF, LB each in 10 ?M concentration)
3. Reverse Transcriptase for Lamp reaction
4. Random primers

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RT-LAMP REV reaction conditions:

REAL-TIME PCR			
Stage	Temperature	Time	Cycles number
Reverse transcription	45 °C	30 min	1
Amplification	60 °C	30 sec	99
		Fluorescence measurement of amplification product after each reaction cycle i real-time	
Product melting	65 ? 95 °C	0,5 °C / 5 sec	
		Fluorescence measurement after each 0.5 °C temperature change.	

RT_LAMP_REV					
					MELTING
	CycleTemp1	CycleTemp1			95°C
		60°C	65°C		
	45°C	odczyt co 30 sek		odczyt co 15 sek	
	30 min	90 cycles			

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Reagent	Volume / 1 rxn	
RT-LAMP REV Isothermal Master Mix	15 ?l	
RT-LAMP REV Primer F3/B3 10 ?M ? 0,5 ?l Primer FIP/BIP 10 ?M ? 2 ?l Primer LF/LB 10 ?M ? 1 ?l	7 ?l	
Super Verte MMLV	1 ?l	
Random primers	1 ?l	
RT-LAMP REV Positive Control / RNA matrix	1 ?l	
Total:	25 ?l	

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Fast RNA purification KIT magnetic particles & 3-Zone (Novazym Polska)


FTA card protocol

FTA SAMPLE PREPARATION:

- 1) Place a piece of FTA card (30-40 mg) in 0.5 mL reagent **3-ZONE**

- 2) Incubate sample for 5 minutes at 15 to 30°C to permit the complete dissociation of nucleoprotein complexes.

- 2) Add 0.2 mL of chloroform per 1 mL of 3-ZONE Reagent.


- 3) Cap sample tubes securely, and shake them vigorously by hand for 15 seconds.


- 4) Incubate tubes at 15 to 30°C for 2 to 3 minutes.

- 5) Centrifuge the samples at 12,000 × g 15 minutes in 2 to 8°C.


- 6) Following centrifugation, the mixture separates into a lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. The volume of the aqueous phase is about 60% of the volume of 3-ZONE Reagent used for homogenization.

RNA PURIFICATION:

- 1) Transfer the aqueous phase to a fresh tube, add 0.5 mL of 70% ethanol and mixed well. Add 25 µl of magnetic beads and incubate 10 min. at room temperature.

- 2) Put the tubes on the magnetic rack and then remove the supernatant.

- 3) Transfer the tubes to an ordinary rack (without a magnet) and wash with 1 mL of 70% EtOH.


- 4) Repeat p.3

- 5) Remove residual ethanol (very carefully). Leave the open tubes for 10-15 minutes in a magnetic rack to evaporate residual ethanol.

- 6) Remove the tubes from the magnetic rack. Resuspend beads in 50 µl of deionised water (free of RNA). Incubate 5 minutes at 75 °C, 800 rpm. Then, open the tubes and incubate for another 5 minutes.

- 7) Transfer the tubes to the magnetic rack. Supernatant, containing RNA, transfer to a new, clean and RNAses free tube.

- 8) Measure of the RNA concentration.