Dane aktualne na dzieñ: 07-05-2024 15:06

Link do produktu: <http://www.novazym.sklep.pl/novabeads-total-rna-purification-kit-100-preparations-3-zone-magnetic-particles-p-437.html>

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

Total RNA purification KIT magnetic particles & 3-Zone (Novazym Polska)

TISSUE PREPARATION:


- 1) Isolate of nucleoprotein complex by using a power homogenizer (BIOPREP-24 or equivalent) in **3-ZONE** (1 mL reagent per each 50-100 mg of tissue). The tissue sample volume should not exceed 10% of the volume of 3-ZONE reagent used for homogenization.


- 2) Incubate the homogenized sample for 5 minutes at 15 to 30°C to permit the complete dissociation of nucleoprotein complexes.

- 2) Add 0.2 mL of chloroform per 1 mL of 3-ZONE Reagent.


- 3) Cap sample tubes securely, and shake them vigorously by hand for 15 seconds.

- 4) Incubate tubes at 15 to 30°C for 2 to 3 minutes.

- 5) Centrifuge the samples at 12,000 × g, 15 minutes in temperature 2 to 8°C.

- 6) Following centrifugation, the mixture separates into a lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. The volume of the aqueous phase is about 60% of the volume of 3-ZONE Reagent used for homogenization.

RNA PURIFICATION:

- 1) Transfer the aqueous phase to a fresh tube, add the same volume of 70% ethanol and mixed well. Add 25 µl of magnetic beads and incubate 10 min. at room temperature.

- 2) Put the tubes on the magnetic rack and then remove the supernatant.

- 3) Transfer the tubes to an ordinary rack (without a magnet) and wash with 1 mL of 70% EtOH.


- 4) Repeat p.3

- 5) Remove residual ethanol (very carefully). Leave the open tubes for 10-15 minutes in a magnetic rack to evaporate residual ethanol.

- 6) Remove the tubes from the magnetic rack. Resuspend beads in 50 µl of deionised water (free of RNA). Incubate 5 minutes at 75 °C, 800 rpm **. Then, open the tubes and incubate for another 5 minutes.

- 7) Transfer the tubes to the magnetic rack. Supernatant, containing RNA, transfer to a new, clean and RNases free tube.

- 8) Measure of the RNA concentration.

* The inhibition during amplification step occurred repeat the p.3 in RNA purification protocol.

** optically. When the elution is carried out by means of dry block without of shaking, the eluates should be

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Fast RNA purification KIT

magnetic particles & 3-Zone (Novazym Polska)


FTA card protocol

FTA SAMPLE PREPARATION:


- 1) Place a piece of FTA card (30-40 mg) in 0.5 mL reagent 3-ZONE

- 2) Incubate sample for 5 minutes at 15 to 30°C to permit the complete dissociation of nucleoprotein complexes.

- 2) Add 0.2 mL of chloroform per 1 mL of 3-ZONE Reagent.


- 3) Cap sample tubes securely, and shake them vigorously by hand for 15 seconds.


- 4) Incubate tubes at 15 to 30°C for 2 to 3 minutes.

shake vigorously

- 5) Centrifuge the samples at 12,000 × g 15 minutes in 2 to 8°C.

- 6) Following centrifugation, the mixture separates into a lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. The volume of the aqueous phase is about 60% of the volume of 3-ZONE Reagent used for homogenization.

RNA PURIFICATION:

- 1) Transfer the aqueous phase to a fresh tube, add 0.5 mL of 70% ethanol and mixed well. Add 25 µl of magnetic beads and incubate 10 min. at room temperature.

- 2) Put the tubes on the magnetic rack and then remove the supernatant.

- 3) Transfer the tubes to an ordinary rack (without a magnet) and wash with mL of 70% EtOH.


- 4) Repeat p.3

- 5) Remove residual ethanol (very carefully). Leave the open tubes for 10-15 minutes in a magnetic rack to evaporate residual ethanol.

- 6) Remove the tubes from the magnetic rack. Resuspend beads in 50 µl of deionised water (free of RNA). Incubate 5 minutes at 75 °C, 800 rpm. Then open the tubes and incubate for another 5 minutes.

- 7) Transfer the tubes to the magnetic rack. Supernatant, containing RNA, transfer to a new, clean and RNases free tube.

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