

Dane aktualne na dzień: 01-05-2025 22:57

Link do produktu: <http://www.novazym.sklep.pl/bacti-cinerator-sterilizer-is800-a-p-3800.html>




## Bacti-cinerator Sterilizer IS800-A

### Opis produktu

# Bacti-cinerator Sterilizer

## IS800-A



### Description:

- ◆ IS800-A Bacti-cinerator Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biosafety cabinets, clean benches. They are also used besides exhaust fans or on mobile vans.

### Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other other severe environment can be easily used.

### Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5 seconds to 7 seconds are needed, because of its highest temperature of 825 °C.
- ◆ Smart, light, easy to handle.
- ◆ Used in anaerobic chamber is allowed.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safer.
- ◆ 12-month warranty.

### Parameters:

| Project                    | Parameters      |
|----------------------------|-----------------|
| Max. Temp. in central zone | 825 ±25         |
| Length of heating zone     | 140mm           |
| Ambient temperature        | 5 ~40           |
| Voltage                    | AC220V~ 50-60Hz |
| Dimension(mm)              | 120×155×180     |

| Project                               | Parameters   |
|---------------------------------------|--------------|
| Max. Out diameter of drug elimination | Φ15mm        |
| Heating time                          | 15min        |
| Fuse                                  | 250V 3A 5×20 |
| Power                                 | 150W         |
| Weight(kg)                            | 1.0          |

Novazym  
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies  
ul. Umultowska 89 C build. A/38, 61-614 Poznań