

Dane aktualne na dzień: 02-05-2025 01:23

Link do produktu: <http://www.novazym.sklep.pl/hybridization-instrument-sh1000-p-4014.html>



Hybridization Instrument SH1000

Opis produktu

Hybridization Instrument

SH1000



Description:

- ◆ SH1000 Hybridization Instrument allows setting constant temperature and fixed time. Compared with traditional manual operation, SH1000 gets a high accuracy about the experiment. It improves people's work efficacy by operating hybridization and denaturizing simultaneously. Temperature is the most part in the experiment of in-situ hybridization and FISH. The instrument may finish the orders that hybridization and denaturizing automatically, according to experimenter's set programs. It assures the research comes along that stick temperature-controlling, appropriate humidity environment, lower cost and shorter time are provided.

Features:

- ◆ Decrease operation time and the hurt from reagent due to FISH processing step.
- ◆ Keep the uniformity of temperature without using 12 slides simultaneously.
- ◆ The design of slide card slot.
- ◆ Cut down on reagent, manpower, capital.
- ◆ Easy to clean the panel.
- ◆ Numeric keyboard is convenient.
- ◆ 45 custom systems.
- ◆ 3 operation models.
- ◆ Internal power supply provided is UL, CSA and CE approved. 12-month warranty.

Parameters:

Project	Parameters
Temp. Range	RT.+5 ~100C
Temp. Accuracy	±0.1C
Heating time	≤3min(from 37C to 95C)
Humidity control	Yes
Voltage	AC 220V/110V(optional)
Weight(kg)	4.5

Project	Parameters
Time range	1min ~ 99h59min
Display Accuracy	0.1C
Cooling time	≤7min(from 95C to 45C)
Power	330W
Dimension(mm)	440×220×120
---	---

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań