

Link do produktu: <http://www.novazym.sklep.pl/manual-headspace-sampler-p-514.html>



Manual Headspace Sampler

Dostępność

Dostępny

Opis produktu

Manual Headspace Sampler

Model: HS-600

Headspace analysis has generally been limited to costly automated systems — until now! This manual sampler enables you to utilise the headspace technique in a quantitative, economical and precise manner. Analyse up to six samples at a time using all popular vial sizes, from 2 - 27ml.

The accurate oven reaches stabilization in around 8 minutes and can heat samples up to 120°C ($\pm 0.5^\circ\text{C}$). Safety is guaranteed with an automatic cut-out that limits the oven temperature to 100°C

Features

- 1? Volatiles in pharmaceuticals
- 2? Flavours analysis
- 3? Alcohol and other toxic compounds in blood
- 4? Screening of volatiles in all type of environmental samples e.g. Soils, waters, plastics, polymers, etc

Specification

- 1? Heating temperature range: RT+5?~100?
- 2? Variable injection: up to 2.5ml
- 3? Temperature accuracy: $\pm 0.5^\circ\text{C}$
- 4? Holds up to 6 vials of: 2, 10, 20ml
- 5? Stabilization time from 25? to 70? with 1ml syringe and 6 empty 20ml vials: 8 minutes

