

Link do produktu: <http://www.novazym.sklep.pl/molecular-diagnostic-kit-hsv-12-hv004b-p-4074.html>



Molecular Diagnostic Kit - HSV 1&2 (HV004B)

Opis produktu

Kit for the simultaneous identification and quantification of Herpes simplex Virus 1 genome by amplification of a HSV1 ? specific sequence within US4 gene and Herpes simplex Virus 2 genome by amplification of a HSV2 ? specific sequence within US4 gene.

Tests: 50

Technique: Real-time PCR