

Link do produktu: <http://www.novazym.sklep.pl/malachite-green-crystal-violet-residue-rapid-test-kit-40-tests-p-4357.html>



Malachite Green, Crystal Violet Residue Rapid Test Kit - 40 tests

Numer katalogowy

RTAN032

Opis produktu

Malachite Green, Crystal Violet Residue Rapid Test Kit Novazym Diagnostics

Target sample: seafood, fish, shrimps

Sensitivity: 2ppb

Tests/kit: 40

Catalogue No.: RTAN032