Dane aktualne na dzień: 01-05-2025 20:38

Link do produktu: http://www.novazym.sklep.pl/vetpcr-c-perfringens-96d-p-623.html



# VetPCR? C. perfringens (96D)

Numer katalogowy

VET-A006-96D

## Opis produktu

Clostridium perfringens PCR Detection Kit

Clostridial necrotic enteritis is a relatively common disease in avian species. Clostridium perfringens is a large Gram negative, sporeforming, anaerobic bacterium that causes necrotic enteritis and subclinical disease in chickens and turkeys. Infection in commercial poultry has been increasing at an alarming rate. Disease is generally caused by type-C strains of C. perfringens producing toxin in the small intestines of birds, resulting in rapped loss of condition and weight loss, lethargic behavior, decreased appetite, and blood stained or undigested food. The toxin and its effects may remain in the system for extended periods of time even after the original bacterial infection has been treated.

VetPCR? C. perfringens Detection Kit is the direct detection of Clostridium perfringens on de basis of a g enetic data base, so it candiagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

#### Characteristics

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20?C. Time-saving and cost-effective.

#### Contents

### Interpretation of the Test Result

Expected PCR product size : 375bp

#### Fig. 1 Result:

Lane M: Brig? Molecular Weight Marker Lane  $1{\sim}2$ : C. perfringens Positive samples

Lane I.C.: Internal control
Lane P: Positive control
Lane N: Negative control