

Novazym http://www.novazym.sklep.pl info@novazym.pl

Dane aktualne na dzień: 02-05-2025 23:06

Link do produktu: http://www.novazym.sklep.pl/vetpcr-lyme-disease-48d-p-754.html



VetPCR? Lyme disease (48D)

Numer katalogowy

VET-B027-48D

Opis produktu

Borrelia burgdorferi PCR Detection Kit

Lyme borreliosis, the most common vector-borne disease in the northern hemisphere, is caused by bacteria Borrelia burgdorferi. Thedisease is multisystemic, affecting mainly the skin, nervous system, heart and joints. In Europe, the vector of the disease is the tick Ixodes ricinus, whereas in the United States of America, two primary tick vectors exist, namely: I. scapularis, and I. pacificus. A number of domestic animals (dogs, horses, sheep and cattle) have contracted Lyme disease. Wild animals (mice, voles, rabbit, deer, racoons, and birds) are known to be reservoir hosts.

VetPCR? Lyme disease Detection Kit is the direct detection of Borrelia burgdorferi on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate, reliable technique.

Characteristics

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20?C. Time-saving and cost-effective.

Contents

Kit	Quantity 48	Quantity 96	Package
VetPCR™ Lyme disease Premixture	1	1	Vial
PCR Internal Control	1	1	Vial
DNase/RNase free water	1	1	Vial
Lyme disease PCR Positive control	1	1	Vial
PCR Negative control	1	1	Vial
Mineral Oil Solution	1	1	Vial
Brig [™] Molecular Weight marker	1	1	Vial
DNA purification kit	50	100	Test

Interpretation of the Test Result

Expected PCR product size : 341bp








Fig. 1 Result:

- Lane M: Brig? Molecular Weight Marker
- Lane 1~2: Lyme disease Positive samples
- Lane I.C.: Internal control
- Lane P: Positive control
- Lane N: Negative control