

Novazym http://www.novazym.sklep.pl info@novazym.pl

Dane aktualne na dzień: 20-05-2024 13:47

Link do produktu: http://www.novazym.sklep.pl/vetpcr-amarginale-96d-p-703.html



# VetPCR? A.marginale (96D)

Numer katalogowy

VET-B001-96D

## Opis produktu

# Anaplasma marginale PCR Detection Kit

Anaplasmosis is caused by a parasite called, Anaplasma marginale. It invades and multiplies in red blood cells. Ticks are biological vectors of A. marginale and spread of the disease is most likely to occur through infected tick bites. Animals sick with Anaplasmosis may have a fever, and be off feed, depressed, dehydrated, show rapid or difficult breathing and have anemia.

VetPCR?A.marginale Detection Kit is the direct detection of Anaplasma marginale on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

#### **Characteristics**

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20?C. Time-saving and cost-effective.

#### Contents

| Kit                                       | Quantity<br>48 | Quantity<br>96 | Package |
|---|----------------|----------------|---------|
| VetPCR™ A.marginale Premixture            | 1              | 1              | Vial    |
| PCR Internal Control                      | 1              | 1              | Vial    |
| DNase/RNase free water                    | 1              | 1              | Vial    |
| A.marginale PCR Positive control          | 1              | 1              | Vial    |
| PCR Negative control                      | 1              | 1              | Vial    |
| Mineral Oil Solution                      | 1              | 1              | Vial    |
| Brig <sup>™</sup> Molecular Weight marker | 1              | 1              | Vial    |
| DNA purification kit                      | 50             | 100            | Test    |

### Interpretation of the Test Result

Expected PCR product size : 117bp



Novazym http://www.novazym.sklep.pl info@novazym.pl



## Fig. 1 Result:

Lane M: Brig? Molecular Weight Marker Lane 1~2: A.marginale Positive samples Lane I.C.: Internal control Lane P: Positive control

Lane N: Negative control