

Dane aktualne na dzień: 05-06-2024 13:30

Link do produktu: <http://www.novazym.sklep.pl/aflatoxin-residue-rapid-test-kit-40-tests-p-4363.html>



Aflatoxin Residue Rapid Test Kit - 40 tests

Numer katalogowy

RTAN058

Opis produktu

Aflatoxin Residue Rapid Test Kit

Novazym Diagnostics

Target sample: edible oil
Sensitivity: 5ppb
Tests/kit: 40
Catalogue No.: RTAN058