

Novazym http://www.novazym.sklep.pl info@novazym.pl

Dane aktualne na dzień: 29-05-2024 06:14

Link do produktu: http://www.novazym.sklep.pl/vetpcr-c-perfringens-96d-p-623.html



VetPCR? C. perfringens (96D)

Numer katalogowy

VET-A006-96D

Opis produktu

Clostridium perfringens PCR Detection Kit

Clostridial necrotic enteritis is a relatively common disease in avian species. Clostridium perfringens is a large Gram negative, sporeforming, anaerobic bacterium that causes necrotic enteritis and subclinical d isease in chickens and turkeys. Infection in commercial poultry has been increasing at an alarming rate. D isease is generally caused by type-C strains of C. perfringens producing toxin in the small intestines of birds, resulting in rapped loss of condition and weight loss, lethargic behavior, decreased appetite, and blood stained or undigested food. The toxin and its effects may remain in the system for extended periods of time even after the original bacterial infection has been treated.

VetPCR? C. perfringens Detection Kit is the direct detection of Clostridium perfringens on de basis of a g enetic data base, so it candiagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 3 hours for detection. Therefore, it is a very fast, accurate and reliable technique.

Characteristics

Ready to use : only DNA template and D.W. are needed. Easy and speed protocol. Stable for 1 year at -20?C. Time-saving and cost-effective.

Contents

Interpretation of the Test Result Expected PCR product size : 375bp

Fig. 1 Result:

Lane M: Brig? Molecular Weight Marker Lane 1~2: C. perfringens Positive samples Lane I.C.: Internal control Lane P: Positive control Lane N: Negative control