

Link do produktu: <http://www.novazym.sklep.pl/agile-plate-illuminator-p-320.html>



Agile Plate Illuminator

Dostępność	Na zamówienie
Czas wysyłki	14 dni
Numer katalogowy	APL-96

Opis produktu

Agile Plate Illuminator

Product specification:

Agile Illuminator provides a simple way to track the loading of reactants by illuminating the target wells which can be single well or a complete/partial row or column of wells, or all 96 wells.

Key features:

- Compatible for 96-well PCR Plate and ELISA Plate
- 3 selectable light intensities
- Applicable with all standard single channel, 8-channel, 16-channel pipettors
- Intuitive and easy instrument with pre-programmed loading patterns
- Water resistant, Compact and Portable

Art. Nr	APL-96
Size	18.8 (L) x 14 (W) x 1.2 (H) cm
Weight	280 g
Vessel Types	96-well PCR Plate and ELISA Plate
Programs	12 preset and user-selectable loading patterns
Light Source	White LED or Red LED (optional)
Operating Voltage	3.5 . 5.3V
Operating Power Consumption	400 mA
Power Supply	100 - 240 V input power, 5 VDC, 1A
Package	Agile Illuminator Unit, Foot pedal and mini-USB