

Link do produktu: <http://www.novazym.sklep.pl/biological-indicator-incubator-minit-3-p-519.html>



Biological Indicator Incubator MiniT-3

Dostępność

Dostępny

Opis produktu

Biological Indicator Incubator Model? MiniT-3

The MiniT-3 Biological Indictaor Incubator is designed for single temperature incubation of the EO Biological Indicators. The thermometer is included with the incubator, provides easy verification of incubation temperature.

Features

1. LCD display, with 24V power for ease of portability
2. Operating Temperature Range: From 5°C above ambient to 100°C
3. Time Range: 1min ~ 99h59min
4. Temp.Control Accuracy: $\leq \pm 0.5^{\circ}\text{C}$
5. Display Accuracy: 0.1°C
6. Heating Time: < 6min (from 20°C to 57°C)
7. Temp.Control Accuracy @57°C: $\leq \pm 0.35^{\circ}\text{C}$
8. Dimensions: 120(w) x 152(D) x 112(H) mm

Specification

1. The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made. User can set the incubating temperature.
2. Matereials: Aluminium block
3. Well Diameter: According to custom-making

Type	Description	Remark
A-MiniT-1	15 x ϕ 10mm	
B-MiniT-1	12 x ϕ 12mm	
C-MiniT-1	6 x ϕ 15mm	
D-MiniT-1	6 x ϕ 16mm	
E-MiniT-1	18 x ϕ 8.8mm	