

Link do produktu: <http://www.novazym.sklep.pl/blok-grzejny-k-30-heating-block-p-127.html>



Blok grzejny K 30 (Heating block)

Dostępność

Na zamówienie

Numer katalogowy

BGK30

Opis produktu

Block for K30 Dry Bath

Type	Parameter	Capacity	Dimension (mm)
A-30	Centrifuge tubes	96×0.2ml	∅108 x 32
B-30	Centrifuge tubes	29×0.5ml and 32×1.5ml	∅108 x 32
C-30	Centrifuge tubes	61×0.5 ml	∅108 x 32
D-30	Centrifuge tubes	41×1.5ml	∅108 x 32
E-30	Centrifuge tubes	41×2.0ml	∅108 x 32
F-30	∅6mm	61	∅108 x 32
H-30	∅12mm	36	∅108 x 50
I-30	∅13mm	36	∅108 x 50
J-30	∅15mm	21	∅108 x 50
K-30	∅16mm	21	∅108 x 50
L-30	∅18mm	17	∅108 x 50