

Dane aktualne na dzień: 01-05-2025 22:29

Link do produktu: <http://www.novazym.sklep.pl/centrifuge-for-microplates-mpc2800-p-4011.html>



Centrifuge for Microplates MPC2800

Opis produktu

Centrifuge For Microplate

MPC2800



Description:

- ◆ MPC2800 is designed for quick spins of samples in PCR plates. Primarily designed for use with 96 or 384 well plates and small volumes, it will accommodate skirted, non skirted and semi skirted styles.

Features:

- ◆ Easy to use, real-time display of all operating information.
- ◆ Countdown mode support
- ◆ DC brushless motor, low noise, small interference, maintenance-free.
- ◆ When opening the cap will automatically stop running.
- ◆ Fast acceleration from standstill to maximum speed fully in just 10 sec.
- ◆ Braking time is short, from the highest speed to a complete standstill in just 4 sec.

Parameters:

Project	Parameters
Rotating speed	2200~2800Rpm
Timing range	15sec ~1min59sec
Max(RCF)	550g
Acceleration time	≤10sec(0 to 2800Rpm)
Brake time	≤4sec(2800 to 0Rpm)
Dimensions	260x240x226mm
Weight	3.3Kg

Project	Parameters
Lid Open Auto Stop	Yes
Short Run	Yes
Capacity	2 PCR plates
Power	45W
Voltage	AC110~220V/50-60Hz
Fuses	250V/3A Φ5×20
---	---