



ChampGel 5000+

Numer katalogowy

CG5000+

Opis produktu

ChampGel 5000+

Specyfikacja:

Cechy indywidualne danego modelu:

System	ChampGel 5000+
Resolution	1392 x 1040 (1.4 megapixels)
Sensor Size (inch)	1/2
Pixel Size	4.65 x 4.65

Wspólne dla wszystkich wersji systemu ChampGel

Pixel Depth	
Lens	Motorized could
Filters	
Trans-White	20 x 30 lighting//wh
Trans-UV	30 1. 254, 302, 2. 21 x 2 3. dimma
Trans-Blue	12 x 12 c
Epi-White	
Epi-UV	

Epi-RGB optional	P.S. Machine
Viewing Window	
Gel operation window	
Software	Lane 1
PC	Not
Screen	Not