

Dane aktualne na dzień: 01-05-2025 23:10

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-dh100-2-p-4018.html>



## Dry Bath Incubator DH100-2

### Opis produktu

# Dry Bath Incubator

DH100-2



## Description:

- ◆ DH100 series high temperature dry bath incubator controlled by microcomputer, use of high purity aluminum material as a heat conduction medium instead of the conditional water bath device. which is easy to use, wide temperature range, high precision and etc. Can be widely used in the sample preservation, reaction of DNA amplification, electrophoresis, degeneration, solidification and other kinds of serum biochemical sample thermostatic incubation process.

## Features:

- ◆ Real-time temperature display, the countdown display.
- ◆ Convenient module replacement, easy to clean and disinfect.
- ◆ Automatic fault detection and alarm function.
- ◆ Built-in over temperature protection device.
- ◆ Temperature deviation calibration function.

## Parameters:

| Project          | Parameters          |
|------------------|---------------------|
| Temp.Range       | RT.+5C -160         |
| Temp uniformity  | ±0.5                |
| Display accuracy | 0.1                 |
| Heating time     | ≤15 min (25 to160 ) |
| Max power        | 500W                |
| Size(mm)         | 260×220×100         |

| Project        | Parameters   |
|----------------|--------------|
| Temp accuracy  | ±0.5 (@40 )  |
| Temp accuracy  | ±1 (@120 )   |
| Time range     | 1-99h59min/∞ |
| Block quantity | 2            |
| Voltage        | 220V 50/60Hz |
| Weight(kg)     | 5.5          |

## Optional Blocks:

| Type | Diameter of test tube | Number of test tube |
|------|-----------------------|---------------------|
| MD01 | 6mm                   | 48                  |
| MD02 | 7mm                   | 48                  |
| MD03 | 10mm                  | 24                  |
| MD04 | 12mm                  | 24                  |
| MD05 | 13mm                  | 24                  |
| MD06 | 15mm                  | 16                  |
| MD07 | 16mm                  | 16                  |
| MD08 | 19mm                  | 12                  |
| MD09 | 20mm                  | 12                  |

| Type | Diameter of test tube  | Number of test tube |
|------|------------------------|---------------------|
| MD10 | 26mm                   | 8                   |
| MD11 | 28mm                   | 4                   |
| MD12 | 40mm                   | 3                   |
| MD13 | 0.2ml centrifuge tube  | 48                  |
| MD14 | 0.5ml centrifuge tube  | 48                  |
| MD15 | 1.5ml centrifuge tube  | 24                  |
| MD16 | 2.0ml centrifuge tube  | 24                  |
| MD17 | 0.2ml well micro plate | 96                  |
| ---  | ---                    | ---                 |

**Novazym**  
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies  
ul. Umultowska 89 C build. A/38, 61-614 Poznań