

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-minit-100h-p-310.html>



Dry Bath Incubator MiniT-100H

Dostępność

Na zamówienie

Numer katalogowy

MINIT100H

Opis produktu

MiniT Series Dry Bath Incubator

Specification

The MiniT dry bath incubator combine a microprocessor-controlled, a fast heating micro-tube block. The Mini series are suitable for a range of uses in life science, molecular biology, environmental and industrial laboratories and a variety of applications including incubation, enzyme reaction, immune-assays. The small form factor allows these to find a place on even the most crowded laboratory benchtop and easily to move.

MODEL	MiniT	MiniT-100	MiniT-100H
Temperature Control Range	RT +5 °C - 80 °C	RT +5 °C - 100 °C	RT +5 °C - 100 °C
Timing Range	0 - 999 min or 0 - 999 sec	0 - 999 min or 0 - 999 sec	1 min - 99 h 59 min
Temperature Control Accuracy	± 0.5 °C	± 0.5 °C	± 0.5 °C
Display Accuracy	0.1 °C	0.1 °C	0.1 °C
Heating Time (20 °C to Max. Temp)	≤ 12 min	≤ 15 min	≤ 15 min
Temp. Control Accuracy @40 °C	0.3 °C	0.3 °C	0.3 °C
Heating Lid	-----	-----	Lid (*)
Heating Lid Temperature Range	-----	-----	RT +5 °C - 100 °C
Dimensions (L x W x H) mm	150 x 110 x 80	152 x 120 x 112	152 x 120 x 115
Net Weight	0.5 kg	0.85 kg	0.9 kg