

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-minit-c-p-462.html>



# Dry Bath Incubator MiniT-C

Dostępność

**Dostępny**

## Opis produktu

The MiniT is a microprocessor-controlled product using advanced thermoelectric technique. Each unit includes a clear plastic cover that provides enhanced temperature uniformity and ramping time. A convenient block lifter is also included, allowing block to be removed from the instruments before they have cooled.

### Features

1. Mini-dimensions: **110 x 150 x 80 mm (Only have the adult palm's sizes)**
2. LCD display. It is easy to setup and use
3. Interchangeable aluminum alloy sample blocks provide versatility and allow easy cleaning and disinfecting
4. Beep-signal / Stop after program completion
5. With 12V power for total portability

### Specification

1. Temperature Control Range: from 5°C above ambient to 80?
2. Timing Range: 1min~99h59min
3. Temperature Control Accuracy:  $\leq \pm 0.3?$
4. Display Accuracy: 0.1?
5. Heating Time:  $\leq 10$  min (from 20? to 80?)
6. Heating Parts: Heater
10. Standard Blocks: See " [Dry Bath Blocks for MiniT](#)"

### Blocks for Dry Bath Incubator MiniT

Type	Dimension	Description	Remark
A-MiniT	0.2ml tube	40×0.2ml	8*0.2ml strip
B-MiniT	0.5ml tube	32×0.5ml	
C-MiniT	1.5ml tube	15×1.5ml	
D-MiniT	2.0ml tube	15×2.0ml	
E-MiniT	cuvette	8-12.5×12.5mm	
F-MiniT	15ml tube	4×15ml	
G-MiniT	50ml tube	2×50ml	