

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-mk200-1a-p-311.html>



Dry Bath Incubator MK200-1A

Dostępność

Na zamówienie

Numer katalogowy

MK200-1A

Opis produktu

Dry Bath Incubator MK200-1A

Specification

The MK200-1A Dry Bath Incubator is a microprocessor-controlled product . With this technique, more stable and accurate temperature control is achieved.

Temperature control range	from 5 °C above ambient to 150 °C
Time required for temperature to increase from 40 °C to 150 °C	≤30 min
Temperature control discrepancy	≤±0.5 °C
Display discrepancy	±0.1 °C
Homogeneity of block temperature	≤± 0.5 °C
Maximum heating power	200W
Dimension	280 x 240 x 110 mm
Time Range	99 h 59 min
Standard Blocks	See " Dry Bath Blocks for MK200-1A"

Key features

- LED display
- Current temperature and count down time display
- Metal blocks free of sample counter-pollution
- Easy replacement, cleaning and disinfecting of metal blocks
- Built-in over-temperature protection
- Calibrating temperature discrepancy
- Buzzer alarm when temperature control ends
- Only heating design and Custom blocks are available