

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-mk2000-2hl-p-297.html>



Dry Bath Incubator MK2000-2HL

Dostępność

Na zamówienie

Numer katalogowy

MK2000-2HL

Opis produktu

Dry Bath Incubator MK2000-2HL

The MK2000-2HL dry bath incubator holds the thermal cycle's heating lid technology, two blocks and touch screen product with a patent. It's the most luxury dry block heating for all over the world.

Specification

Temperature Control Range	RT+5 °C ~ 100 °C
Stability @ 100 °C	±0.4 °C
Uniformity	
within the block @ 37 °C	±0.2 °C
across similar blocks @ 37 °C	±0.3 °C
Heating time (25 °C to 100 °C)	≤15min
Cooling time (100 °C to 40 °C)	≤20min
Display Accuracy / Display	0.1 °C / Touch Screen
Timing Range	1sec ~ 99h59min
Block Quantity	2
Dimensions L x W x H (mm)	380 x 210 x 160
Weight (kg)	5.5
Standard Blocks	See "Blocks for MK2000-2HL Dry Bath Incubator"

Key features

- Exceptional temperature uniformity and accuracy -- due to advanced temperature control combined with high quality, precision engineered block providing excellent thermal contact
- Hold a heating lid, 60 user programmable settings
- Touch screen allows for easy to read and easy programming and easy to use. It can display the real time temperature
- Programmable system and program link functions, it can be run for some simple PCR tests
- Wide range of interchangeable aluminium blocks
- Block extraction tool is supplied, allowing blocks to be removed easily
- A stainless plate can press on the tubes lid and avoid the tubes lid opening when user incubate at a high temperature