

Link do produktu: <http://www.novazym.sklep.pl/dry-block-heating-for-testing-milk-antibiotic-minis-p-515.html>



Dry Block Heating for Testing Milk Antibiotic - MiniS

Dostępność

Dostępny

Opis produktu

Dry Block Heating for Testing Milk Antibiotic

Model: MiniS

The MiniS dry block heating is designed for testing milk for antibacterial substances such as antibiotics and sulphonamide residues. It includes a digital thermometer and timer. Heater is either AC or DC compatible. AC adapter is included, DC adapter for vehicle must be purchased separately.

Features

1. A personal dry block heating and for digital dry block heating
2. It pretty much fits anywhere on your lab bench
3. It has a very simple touch pad digital display for temperature and time selection
4. It's only 4-inch wide
5. High quality -- low price

Code

Description

AS-01250-00	MiniS DC12V 20W
AS-01251-01	A-MiniS block 3--8 x 0.2ml
AS-01251-02	B-MiniS block 2--8 x 0.2ml + 5 x Ø11mm tubes

Specification

- | | |
|--|-----------------|
| 1. Temperature Control Range: | RT+5? ~ 80? |
| 2. Stability @ 80?: | ±0.5? |
| 3. Uniformity: within the block @ 40?: | ±0.3? |
| 4. Heating time (25? to 80?): | ≤12min |
| 5. Display Accuracy / Display: | 0.1? / LCD |
| 6. Timing Range: | 1min ~ 99h59min |
| 8. Dimensions L x W x H (mm): | 135 x 99 x 33.5 |