

Dane aktualne na dzień: 01-05-2025 22:54

Link do produktu: <http://www.novazym.sklep.pl/e-cooler-model-ec-01-p-169.html>



E-Cooler, model EC-01

Dostępność

Dostępny

Numer katalogowy

EC-01

Opis produktu

E-COOLER
EC-01

[Specification](#)



Key features

- Environment Friendly: Ice-free, just plug-in and press the power button; save energy more than 80%, every eCooler diminishes CO2 emission 2.67MT/year
- Portable: Cooling box keeps 4°C below for 25min without cooler
- Ultra rapid Cooling: 10 min from RT to below 4°C, and maintains 0 ± 0.5°C
- Organized: Easy to organize samples
- Dry: Free your hands and samples from wet ice
- Fits: For all standard microcentrifuge tubes
- Quiet: Ultra low noise (19 dB-A)