

Link do produktu: <http://www.novazym.sklep.pl/ecp3000-high-voltage-multi-usage-power-supply-p-343.html>



ECP3000 High-Voltage Multi-Usage Power Supply

Dostępność

Na zamówienie

Numer katalogowy

ECP3000

Opis produktu

Specification

- Size: 30×26×12(cm)
- Weight: 5kg
- Output Range: 20--3000V, 1--200mA, 1--200W.
- Output jacks: Four sets of output jacks are provided to facilitate connection of up to 4 identical electrophoresis cells simultaneously
- Timer control: 1min--99hr, 59min
- Display: Backlit LCD(128×64pixel)
- Control: Intellectual microprocessor control unit

Key features

- Programmable constant voltage, constant current or constant
- Power operation with automatic crossover
- Safety features include automatic detection of no-load, short circuit, rapid resistance change, ground leak and system overheating
- Ten general electrophoresis methods stored function
- Automatic memory
- Standard and timer operation function
- Automatic completion of a run interrupted by power failure
- Intellectual display of Constant-Voltage, Current, Power
- Adjustable LCD display contrast
- Optional functions (Recovery of power failure, Non-proliferation electric field, V-H?run step by step)