

Link do produktu: <http://www.novazym.sklep.pl/ecpt3000-high-voltage-multi-usage-power-supply-p-344.html>



ECPT3000 High-Voltage Multi-Usage Power Supply

Dostępność

Na zamówienie

Numer katalogowy

ECPT3000

Opis produktu

Specification

- Size: 30×26×12(cm)
- Weight: 5kg
- Output Range: 20--3000V, 1--200mA, 1--200W.
- Output jacks: Four sets of output jacks are provided to facilitate connection of up to 4 identical electrophoresis cells simultaneously
- Timer control: 1min--99hr, 59min
- Display: Backlit LCD(128×64pixel)
- Control: Intellectual microprocessor control unit

Key fetures

- Overload, no load and current leakage protection function
- Ten general electrophoresis methods stored function
- Automatic memory
- Gel temperature tested function (need to purchase relevant temperature sensor separately)
- Standard and Timer operation, auto-switching function
- Automatic completion of a run interrupted by power failure
- Intellectual display of Constant-Voltage, Current or Power
- Temperature can be controlled during the electrophoresis gel experiments through the touch key.
- Standard RS232 connector, which is connected with a PC, analyzing the gel data and temperature variation.
- Optional functions (Recovery of power failure, Non-proliferation electric field, V-H?run step by step)