



Numer katalogowy

4831T-D

Thermal Cycler (PCR) / Gene-Explorer Touch

This series product adopts 8 pieces custom-made US Marlow peltier, Max. ramping rate 5 $^{\circ}$ C/s, cycle times $\leq$ 600000. More advanced technologies applied on the product: Windows operating system; color touch-screen; double blocks can run independently; replaceable block; control by PC function; print function; extra large data storage and USB interfaces. All the above functions make it demand high need experiments.

>>>>>>>>>> Gene-Explorer Touch.pdf