

Link do produktu: <http://www.novazym.sklep.pl/gene-explorer-touch-6021t-p-846.html>



Gene-Explorer Touch (6021T)

Numer katalogowy

6021T-B

Opis produktu

Thermal Cycler (PCR) / Gene-Explorer Touch

This series product adopts 8 pieces custom-made US Marlow peltier, Max. ramping rate 5?, cycle times ?600000. More advanced technologies applied on the product: Windows operating system; color touch-screen; double blocks can run independently; replaceable block; control by PC function; print function; extra large data storage and USB interfaces. All the above functions make it demand high need experiments.

>>>>>>>>>[Gene-Explorer Touch.pdf](#)