

Link do produktu: <http://www.novazym.sklep.pl/genie-centrifuge-p-481.html>



Genie Centrifuge

Dostępność

Dostępny

Numer katalogowy

OP-FUGE

Opis produktu

Reliable and easy-to-use, the Prism? Mini is compact and economical, an excellent accessory for the Genie? II. The Prism? Mini can be used for a wide range of molecular biology separations and quick spins. The centrifuge was designed to be extremely compact with an exceptionally small footprint. It includes 2 quick release interchangeable rotors (8 place 1.5/2.0ml rotor and 4 place PCR strip rotor), which spin up to 6000 RPM?s. The electronic brake and self opening lid allow for quick deceleration to reduce handling time. The Prism? Mini?s maintenance-free drive is remarkably quiet and the unique design allows for excellent air-flow, to protect temperature-sensitive samples.

- New rotor designs for increased sample capacity
- Self opening lid enables single-handed operation
- Quick release rotor hub allows for easy rotor exchange
- Built in braking system for rapid sample processing

