

Link do produktu: <http://www.novazym.sklep.pl/genie-tubes-small-pack-size-50-strips-of-8-tubes-p-478.html>



Genie tubes Small pack size - 50 strips of 8 tubes

Dostępność

Dostępny

Numer katalogowy

OP-0008-50

Opis produktu

Genie[®] II uses a proprietary tube strip that maximises optical and thermal efficiencies with a locking cap providing a closed-tube system. The strips have the following advantages:

- Seal-and-lock mechanism to prevent contamination
- Individually capped
- Non-fluorescent and optically clear
- Wings for ease of handling
- Raised rim for foil sealing
- Each strip has 8-tubes with a working volume of 20 ? 150 ?l

BRAK
ZDJĘCIA



BRĄK
ZDJĘCIA



(strips in holder)