

Link do produktu: <http://www.novazym.sklep.pl/glass-bead-sterilizer-gbs-3000a-p-299.html>



Glass Bead Sterilizer GBS-3000A

Dostępność

Na zamówienie

Numer katalogowy

GBS-3000A

Opis produktu

Glass Bead Sterilizers GBS-3000A

The Glass Bead Sterilizers are a Quick, Easy and Accurate alternative to traditional methods of sterilization. Once switched 'ON', the units are ready to use within 15 minutes and ensure total sterilization by destruction of all Micro - Organisms within 10 seconds.

Glass Bead Sterilizers are incorporated with an imported thermostat, which maintains the temperature with 240 °C - 260 °C. These units, having a stainless steel body, are compact enough to be placed on any Laminar Air Flow Work benches or any other worksite tables in a clean room atmosphere.

Specification

Model	GBS-3000A	GBS-3000B
Crucible size	40 mm Φ80 mm Depth	40 mm Φ140 mm Depth
Preset Temperature Range	240 °C – 260 °C	240 °C – 260 °C
Initial Stabilization Time	10 minutes	15 minutes
Weight of Glass Bead	40g	150g
Material of Lid	Stainless	Stainless
Wattage	50 W Single Phase	60 W single phase
Input	220 V / 50 Hz or 110 V / 50 Hz	220 V / 50 Hz & 110 V / 50 Hz
NOTE: Temperature Range can be set according to customers requirements		