

Dane aktualne na dzień: 01-05-2025 22:59


Link do produktu: <http://www.novazym.sklep.pl/hematocrit-centrifuge-hwlab-p-3810.html>



Hematocrit centrifuge HWLAB

Opis produktu

Hematocrit centrifuge HWLAB



Description

HWLAB® Hematocrit centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 500r/min to 12000r/min
- Timer, from 15s to 99min, or continuous operation
- For 24x75mm hematocrit tubes
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time

Technical specifications

Description	HW24H
Maximum rcf	13680xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	24 x 75mm hematocrit tubes
Speed	From 500r/min to 12000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	212x270x158mm
Weight	4.5kg



Ordering information: HWLAB® Clinical centrifuge

Model	Cat.NO	Destription
 HW24H	C-11000	Widely ues in clinical labs and donation centers
 HW24H.1	C-11001	For 24x75mm hematocrit tubes
 HW24H.2.1	C-11002	75mm hematocrit tubes, non-heparinized
 HW24H.2.2	C-11003	75mm hematocrit tubes, sodium-heparinized