

Dane aktualne na dzień: 04-05-2024 05:10

Link do produktu: <http://www.novazym.sklep.pl/incubator-for-microplates-dh400-p-892.html>



Incubator for Microplates DH400

Opis produktu

Incubator For Microplates

DH400



Description:

- ◆ DH400 incubator for microplates is a kind of microprocessor controlled and PID controlled incubator. It is compact, versatile, reliable and easy to use. Mainly used in elisa plates(96/384 wells) or Tissue culture plates(24/48/96 wells,etc) for samples cultivation of incubation .

Application:

Cytochemistry---for in situ reactions.
Immunochemistry---for immunofermentative reactions.
Biochemistry---for enzyme and protein analysis.
Molecular biology---for micro array analysis.

Features:

- ◆ LCD displays.
- ◆ Easy to operate with one-touch knob.
- ◆ Microprocessor and temperature controlled with good liner, small fluctuations.
- ◆ Double heating system ensures the stability of incubation.
- ◆ Internal power supply provided is UL, CSA and CE approved.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Temp.Range	Ambient +5 ~80
Temp. Accuracy	±0.5
Temp. Uniformity	±0.5
Max. height of plate	22 mm
Voltage	AC 110V~240V
Dimensions(mm)	280×270×140

Project	Parameters
Time range	1min ~ 99h59min
Display Accuracy	0.1
Heating time	≤20min (from 25 to 80)
Capacity	2 plate
Power	120W
Weight(kg)	4.0

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań

