

Dane aktualne na dzień: 01-05-2025 22:21

Link do produktu: <http://www.novazym.sklep.pl/incubator-for-microplates-dh500-p-4013.html>



Incubator For Microplates DH500

Opis produktu

Incubator For Microplates

DH500



NEW in 2016

Description:

- ◆ DH500 incubator for microplates which PID fuzzy control technology can accurate to ensure that the temperature control precision and automatically adjust the heating rate, reduce waiting time. It is the ideal automation tool for samples hatching, catalytic reaction process of preserving, blending, save and etc.

Features:

- ◆ Easy to set up and use, all information real-time display and showing set up operation, convenient to observe equipment running status.
- ◆ Supports standard microplates and deep well plates.
- ◆ Automatic preheating function.
- ◆ Automatic power recovery function.
- ◆ Temperature calibration function.
- ◆ Built-in software and hardware over temperature protection device, will use more reliable.

Parameters:

Project	Parameters
Temp setting range	0°C~80°C
Temp control range	RT+5°C~80°C
Time range	1min ~ 99h59min/∞
Temp accuracy	±0.5°C*
Display accuracy	0.1°C
Temp uniformity	±0.5°C*
Heating time	≤10min(from 25°C to 80°C)*

Project	Parameters
Auto preheating	Yes
Auto resume to run	Yes
Capacity	4 Microplates or Deep well plates
Power supply	300W
Voltage	AC220V/50-60Hz
Size	340x320x200 mm
Weight	6.5Kg

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań

