

Dane aktualne na dzień: 01-05-2025 20:18

Link do produktu: <http://www.novazym.sklep.pl/incubator-shaker-for-microplates-ms100-p-895.html>



Incubator-shaker for Microplates MS100

Opis produktu

Shaker For Microplates

MS100



Description:

- ◆ MS100 Shaker for Microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for elisa plates(96/384 wells), tissue culture plates(24/ 48/96 wells).

Features:

- ◆ LCD display system status and parameters.
- ◆ Stable and reliable shaking.
- ◆ Easy to operate with one touch knob.
- ◆ Automatic fault detection and buzzer alarm function.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Time Range	1min ~ 99h59min
Mixing Orbit	3mm
Sample capacity	4 plates
Voltage	AC110V~240V 50/60 Hz
Weight(kg)	4.0

Project	Parameters
Speed Range	200~1600 Rpm
Speed step value	10Rpm
Power	45W
Dimension(mm)	280×270×140
---	---

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań