

Link do produktu: <http://www.novazym.sklep.pl/jy-cx2b-nucleic-acid-sequencing-p-358.html>



JY-CX2B Nucleic Acid Sequencing

Dostępność

Na zamówienie

Numer katalogowy

JY-CX2B

Opis produktu

JY-CX2B Nucleic Acid Sequencing

JY Nucleic Acid sequencing series products are designed for the AFLP analysis, difference display, SSR and SSCP research activities, also, they are fit for heteroduplex, DNA foot printing, RNase protection test, S1 Nucleic acid enzyme figures and Primer Extension experiments. Recommend matched power supply: JY3000/JY5000/ JY3000+

Main specification of JY-CX1B/JY- CX2B /JY- CX3B Vertical Electrophoresis Series

Model	JY-CX1B	JY-CX2B	JY-CX3B
Gel size (W×L)mm	180×420	330×420	330×320
Buffer volume (ml)	500	750	750
Sample volume	33,40	50,62,74	50,62,74

Advantages

- Simple operation, only four fasten knobs needed.
- Special heat dissipation processed aluminum plate, firmly cling to the gel plates, which is effectively guaranteed the gel temperature more even, ensuring the electrophoresis trip more straight.
- The discharge pipe for the upper and bottom buffer, which makes the filling and discharge more convenient and simple.
- 50% buffer saved., compare with the traditional tank.

Standard Accessory

Description	Quantity
Body Tank	1
Notched glass plates	2
Square glass plates	2
Shark comb	2
Edging strip	4
Clips	8
Sealing strip	1
Electrical Cables	1