


Link do produktu: <http://www.novazym.sklep.pl/jy-jx5l-rapid-ssr-vertical-electrophoresis-system-p-356.html>



JY-JX5L Rapid SSR Vertical Electrophoresis System

Dostępność

Na zamówienie

Numer katalogowy

JY-JX5L

Opis produktu

JY-JX5L Rapid SSR Vertical Electrophoresis System

New Rapid SSR Vertical Electrophoresis System is one new generation products, it is unique in China. From the appearance point of view, Rapid SSR Vertical Electrophoresis System keeps the advantages of "Sandwich" core with abundant buffer, avoiding the buffer and gel leakage fundamentally, no need of sealing strip, simple to operate.

New Rapid SSR Vertical Electrophoresis System is suitable for PAGE and protein analysis electrophoresis. Four operation screws ensuring the gel caster more convenient and simple, safety operation design conception, following are the series products of Rapid SSR Vertical Electrophoresis System. Recommend matched power supply: JY300C/JY600C/ JY1000C

Main specification of JY-JX5/ JY-JX7/JY-JX5L Vertical Electrophoresis Series

Model	JY-JX5	JY-JX5L	JY-JX7
Gel size (W×L)mm	186×95	186×205	230×95
Buffer volume (ml)	1800	3500	2000
Sample volume	25,40,52	25,40,52	54

JY-JX5/JX5L Standard Accessory

Description	Quantity
Body Tank	1
Notched glass plates	4
1mm bonded glass spacer	4
1mm thick 25wells comb	2
1mm thick 40 wells comb	2
1mm thick 52 wells comb	2
Syringe	1
Electrical Cables	1