


Link do produktu: <http://www.novazym.sklep.pl/jy-jx7-rapid-ssr-vertical-electrophoresis-system-p-355.html>

JY-JX7 Rapid SSR Vertical Electrophoresis System

Dostępność

Na zamówienie

Numer katalogowy

JY-JX7

Opis produktu

JY-JX7 Rapid SSR Vertical Electrophoresis System

New Rapid SSR Vertical Electrophoresis System is one new generation products, it is unique in China. From the appearance point of view, Rapid SSR Vertical Electrophoresis System keeps the advantages of "Sandwich" core with abundant buffer, avoiding the buffer and gel leakage fundamentally, no need of sealing strip, simple to operate.

New Rapid SSR Vertical Electrophoresis System is suitable for PAGE and protein analysis electrophoresis. Four operation screws ensuring the gel caster more convenient and simple, safety operation design conception, following are the series products of Rapid SSR Vertical Electrophoresis System. Recommend matched power supply: JY300C/JY600C/ JY1000C

Main specification of JY-JX5/ JY-JX7/JY-JX5L Vertical Electrophoresis Series

Model	JY-JX5	JY-JX5L	JY-JX7
Gel size	186×95	186×205	230×95
(W×L)mm			
Buffer volume	1800	3500	2000
(ml)			
Sample volume	25,40,52	25,40,52	54

JY-JX7 Standard Accessory

Description	Quantity
Body Tank	1
Notched glass plates	4
1mm bonded glass spacer	4
1mm thick 54wells comb	4
Syringe	1
Electrical Cables	1