

Link do produktu: <http://www.novazym.sklep.pl/jy-scz6-rapid-ssr-vertical-electrophoresis-system-p-351.html>



JY-SCZ6 Rapid SSR Vertical Electrophoresis System

Dostępność

Na zamówienie

Numer katalogowy

JY-SCZ6

Opis produktu

JY-SCZ6 Rapid SSR Vertical Electrophoresis System

JY series Doubled Vertical Protein Tanks offers users different gel sizes for chooses, considering the different demands for the experiments, this unit utilizes the most advanced design features to provide convenient ease of use with high resolution separations. Ideal for the 2-D electrophoresis. Recommend matched power supply: JY300C/JY600C/JY 1000C

Main specification of JY-SCZ6 / SCZ7/ JY-SCZF Vertical Electrophoresis Series

Model	JY-SCZ6	JY-SCZ7	JY-SCZF
Gel size (W×L)mm	160×180	230×150	306×95
Buffer volume (ml)	600	1800	3500
Sample volume	21,26	54	102

JY-SCZ6 Standard Accessory

Description	Quantity
Body Tank	1
Notched glass plates	4
1mm bonded glass spacer	2
1.5mm bonded glass spacer	2
1mm thick 35 wells comb	2
1.5mm thick 35 wells comb	2
1mm thick 21 wells comb	2
1.5mm thick 21wells comb	2
1mm and 1.5mm thick sealing strip	4
Clips	12
Electrical Cables	1