

Link do produktu: <http://www.novazym.sklep.pl/jy-scz9-medium-protean-electrophoresis-cell-p-350.html>



JY-SCZ9 Medium-Protean Electrophoresis Cell

Dostępność

Na zamówienie

Numer katalogowy

JY-SCZ9

Opis produktu

JY-SCZ9 Medium-Protean Electrophoresis Cell

JY-SCZ8/SCZ9 vertical electrophoresis tank is intended for the Protein purification and analysis, ideal for the dimensional electrophoresis, DGGE and preparative gel electrophoresis. JY-SCZ8/SCZ9 is equipped with special processing technology aluminum plates, effectively guaranteed the heat dissipation more even, keeping a stable PH value during the run process. Non-proliferation electric field influence, guaranteed the electrophoresis strip more clear and straight. Modular double hinged core cassettes, matched with the basing gel caster, ensuring the whole operation process more convenience. Retractable connectors compatible with most major power supplies. For the temperature control, JY-SCZ8/SCZ9 can be connected with a outer constant temperature circulator. Safety operation design, power will be cut off automatically when remove the lid. Recommend matched power supply: JY300C/JY600C/ JY1000C

Main specification for JY-SCZ9

Gel size: 120 × 145 mm
Sample volume: 16, 23
Buffer Volume: 1 200 ml

JY-SCZ9 Standard Accessory

Description	Quantity
Body Tank	1
Notched glass plates	4
1mm bonded glass spacer	2
1.5mm bonded glass spacer	2
1mm thick 16 wells comb	2
1.5mm thick 16 wells comb	2
1mm thick 23 wells comb	2
1.5mm thick 23 wells comb	2
Gel caster	1
Water pipe	3
Electrical Cables	1