

Link do produktu: <http://www.novazym.sklep.pl/jy-sp11-horizontal-electrophoresis-p-207.html>



JY-SP11 Horizontal Electrophoresis

Dostępność

Na zamówienie

Numer katalogowy

JY-SP11

Opis produktu

JY-SP11 Horizontal Electrophoresis

Specification:

JY-SP11 is specially designed to meet the most customer and market requirements, because of the sample quantity demands, our engineers took great effort to make market surveys and received many feedbacks from the terminal user. JY-SP11 is usually served for the currently electrophoresis available, features with the DNA electrophoresis screening. A total of several rows of 96-sample can be done simultaneously.

Key features:

1. Gel Size: 245×250 mm
2. Sample volume: 52
3. Buffer Volume: 1 700 mL
4. Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
5. Unfolded electrodes make the repairing and parts replacement easily and quickly
6. The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking
7. Auto-switching when lid removal

8. Recommend matched Power Supply: JY300C / JY600C

JY-SP11 Standard Accessory:

DESCRIPTION	QUANTITY
Electrophoresis cell	1
Gel tray	1
Gel caster	1
Gel caster slide	1
1.5mm thick 52-well comb	7
Electrical Cables	1