

Link do produktu: <http://www.novazym.sklep.pl/jy-sp3-rapid-electrophoresis-system-p-204.html>



JY-SP3 Rapid Electrophoresis System

Dostępność

Na zamówienie

Numer katalogowy

JY-SP3

Opis produktu

JY-SP3 Rapid Electrophoresis System

Specification:

JY-SP3 Rapid electrophoresis system is intended for quick check of samples, which is served for the PCR check, the whole experiment process can be finished within 15 minutes. Meanwhile, the 96-well samples is matched with the PCR board well, convenient for the PCR requirements. For a JY-SP3, the gel casting, running and analysis are all performed in the same ultra compact unit. Buffer and gel volumes have been kept to a minimum and the parallel electrode arrangement allows ultra efficient current transfer, enabling resolution to be completed within 15 minutes. Black bottom of cell enable the observation more clear when loading the samples and running. Safety design for the electricity cutting off when removes the lid.

Key features:

1. Gel Size: 85×200 mm
2. Sample volume: 12 x 10 well, 12 x 18 well
3. Buffer Volume: 50 mL
4. Fast electrophoresis running speed, the whole process could be finished within 15 mins
5. Ultra compact
6. Low buffer volumes
7. Multichannel pipette compatible combs for speed loading
8. Auto-switching when lid removal

9. Recommend matched Power Supply: JY300C / JY600C

JY-SP3 Standard Accessory:

DESCRIPTION	QUANTITY
Electrophoresis cell	1
Sample comb shelf	1
Gel Barrie	2
Gel Tray	1
Electrical Cables	1