

Link do produktu: <http://www.novazym.sklep.pl/jy-sp7-cellulose-acetate-membrane-electrophoresis-system-p-205.html>



JY-SP7 Cellulose Acetate Membrane Electrophoresis System

Dostępność

Na zamówienie

Numer katalogowy

JY-SP7

Opis produktu

JY-SP7 Cellulose Acetate Membrane Electrophoresis System

Specification:

JY-SP7 is updated from the JY-SP7A, with double-bridges style design, which is increased the samples in two times compare with the JY-SP7A, proof-structure lid, enable the experiment volatile liquids not drop on the Cellulose Acetate Membranes. JY-SP7 equipped with two sets of adjustable beam for controlling the run distance.

Key features:

1. Bridge Size 1 : 66×250 mm x2 (W×L)
2. Bridge Size 2 : 82×250 mm x2 (W×L)
3. Buffer Volume: 1 000 ml
4. Size of Cellulose Acetate Membrane: 70×90 mm
5. Recommend matched Power Supply: JY600C / JY600

JY-SP7 Standard Accessory:

DESCRIPTION	QUANTITY
Electrophoresis cell	1
Adjustable beam	4
Retractable connectors	1

Cellulose Acetate Membrane(hard)	10
Cooling pipe	2
Adjustable slide leg	4
Electrical Cables	1