

Link do produktu: <http://www.novazym.sklep.pl/jy-spat-horizontal-electrophoresis-p-208.html>



JY-SPAT Horizontal Electrophoresis

Dostępność

Na zamówienie

Numer katalogowy

JY-SPAT

Opis produktu

JY-SPAT Horizontal Electrophoresis

Specification:

The JY-SPAT is the smallest unit in the range, designed for low to medium numbers of samples. The small gel size is 48x75 mm, and combs ranging from 3 up to 9 samples.

Simply by altering the combs, this compact unit is capable of resolving up to 9 different samples, prepping 1 mL of sample or separating sample bands over a distance of 3 cm. Injection molded construction durable leak-proof environment for complete safety and long life.

Key features:

1. Gel Size: 48×75 mm
2. Sample volume: 3, 5, 9
3. Buffer Volume: 150 mL
4. Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
5. Unfolded electrodes make the repairing and parts replacement easily and quickly
6. The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking

7. Auto-switching when lid removal

8. Recommend matched Power Supply: JY300C/JY300

JY-SPAT Standard Accessory:

DESCRIPTION	QUANTITY
Electrophoresis cell	1
Gel tray	1
Gel caster	1
3-sample, 1mm thick comb	2
3-sample, 1.5mm thick comb	2
5-sample, 1mm thick comb	2
5-sample, 1.5mm thick comb	2
9-sample, 1mm thick comb	2
9-sample, 1mm thick comb	2
Electrical Cables	1