

Link do produktu: <http://www.novazym.sklep.pl/jy-spbt-horizontal-electrophoresis-p-209.html>



JY-SPBT Horizontal Electrophoresis

Dostępność

Na zamówienie

Numer katalogowy

JY-SPBT

Opis produktu

JY-SPBT Horizontal Electrophoresis

Specification:

The JY-SPBT is for separation and identification of DNA and RNA fragments, PCR products and synthetic oligonucleotides. With external gel caster and gel tray, the gel making process is easily to be finished. Injection mould construction - durable leak-proof environment for complete safety and long life.

Key features:

1. Gel Size: 78×100 mm
2. Sample volume: 10, 15
3. Buffer Volume: 400 mL
4. Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
5. Unfolded electrodes make the repairing and parts replacement easily and quickly
6. The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking
7. Auto-switching when lid removal
8. Recommend matched Power Supply: JY300C / JY300 / JY600C

JY-SPBT Standard Accessory:

DESCRIPTION	QUANTITY
Electrophoresis cell	1
Gel tray	1
Gel caster	1
3-sample, 1mm thick comb	2
3-sample, 1.5mm thick comb	2
5-sample, 1mm thick comb	2
5-sample, 1.5mm thick comb	2
9-sample, 1mm thick comb	2
9-sample, 1mm thick comb	2
Electrical Cables	1