

Link do produktu: <http://www.novazym.sklep.pl/jy-spct-horizontal-electrophoresis-p-210.html>



# JY-SPCT Horizontal Electrophoresis

Dostępność

**Na zamówienie**

Numer katalogowy

**JY-SPCT**

## Opis produktu

### Specification

The JY-SPCT is ideal for restriction fragment analysis, sample prep or checking of high numbers of samples, also it is suitable for the nucleic acid purification and preparation.

The JY-SPCT offers a wide degree of versatility. Four tray options are available - 120×120 mm, 120×60 mm, 60×120 mm and 60×60 mm, from the mini-size to the medium-size, which meet the most people's choice.

### Key features

- Gel Size: 60×60 mm, 60×120 mm, 120×60 mm, 120×120 mm
- Sample volume: 3, 6, 8, 11, 13, 18, 25
- Buffer Volume: 550 mL
- Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
- Unfolded electrodes make the repairing and parts replacement easily and quickly
- The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking
- Auto-switching when lid removal
- Recommend matched Power Supply: JY300C / JY600C

### JY-SPCT Standard Accessory

DESCRIPTION	QUANTITY
Cell with electrodes	1
120×120 mm gel tray	1
120×60 mm gel tray	1
60×120 mm gel tray	1
60×60 mm gel tray	2
Gel caster	1
1.5 mm thick 3 wells comb	1
1 mm thick 6&13 wells comb	2
1 mm thick 8&18 wells comb	2
1 mm thick 11&25 wells comb	2
Electrical Cables	1