

Link do produktu: <http://www.novazym.sklep.pl/jy-spdt-horizontal-electrophoresis-p-211.html>



JY-SPDT Horizontal Electrophoresis

Dostępność

Na zamówienie

Numer katalogowy

JY-SPDT

Opis produktu

Specification

The same as the general horizontal electrophoresis cell, JY-SPDT is ideal for the molecular biology research, especially for the nucleic acid purification and preparation. As to the Cells itself, JY-SPDT adopts with injection molding technology of high-strength polycarbonate. Compared with JY-SPCT, JY-SPDT has relative more bigger than that of JY-SPCT, three kinds of gel size are available for it, because of some special experiment needs, JY-SPDT offers the great convenient for the experiments operator.

Key features

- Gel Size: 140×140 mm, 140×70 mm, 70×70 mm
- Sample volume: 3, 6, 8, 13, 14, 19, 29
- Buffer Volume: 900 mL
- Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
- Unfolded electrodes make the repairing and parts replacement easily and quickly
- The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking
- Auto-switching when lid removal
- Recommend matched Power Supply: JY300C / JY600C

JY-SPDT Standard Accessory

DESCRIPTION	QUANTITY
Electrophoresis cell	1
140×140mm gel tray	1
70×70mm gel tray	2
140×70mm gel tray	1
Gel caster	1
1.5mm thick 3 wells comb	1
1mm thick 6 & 14 wells comb	1
1mm thick 8 & 19 wells comb	2
1mm thick 13 & 29 wells comb	2
Electrical Cables	1
Water saving module	2