

Link do produktu: <http://www.novazym.sklep.pl/jy-spft-horizontal-electrophoresis-p-213.html>



# JY-SPFT Horizontal Electrophoresis

Dostępność

**Na zamówienie**

Numer katalogowy

**JY-SPFT**

## Opis produktu

### Specification

The JY-SPFT is designed for ease of use and safe operation, the whole horizontal electrophoresis systems and accessories are built to withstand the rigours of everyday use.

For leak proof performance, the cell itself and UV transparent gel trays are made of high polycarbonate with injection moulding. Cassettes protect the electrodes and allow for easy replacement. The JY-SPFT is equipped with three kinds of mostly used trays and combs, ideal for the DNA and RNA agarose gel electrophoresis currently available; also, it is suitable for the general electrophoresis experiments.

### Key features

- Gel Size: 130×130 mm, 100×78 mm, 48×75 mm
- Sample volume: 3, 5, 10, 13, 15, 26
- Buffer Volume: 600 mL
- Abundant buffer not only guarantee the cooling effects, but also keep the pH Value stable during the whole experiment process
- Unfolded electrodes make the repairing and parts replacement easily and quickly
- The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids the electricity leaking
- Auto-switching when lid removal

Recommend matched Power Supply: JY300C / JY600C

### JY-SPFT Standard Accessory

| DESCRIPTION               | QUANTITY |
|---------------------------|----------|
| Cell with electrodes      | 1        |
| 130×130mm gel tray        | 1        |
| 78×100mm gel tray         | 1        |
| 48×75mm gel tray          | 2        |
| Gel caster                | 1        |
| 1mm thick 13 wells comb   | 2        |
| 1.5mm thick 13 wells comb | 2        |
| 1mm thick 26 wells comb   | 2        |
| 1.5mm thick 26wells comb  | 4        |
| 1mm thick 10 wells comb   | 2        |
| 1mm thick 15 wells comb   | 2        |
| 1.5mm thick 3 wells comb  | 2        |
| 1mm thick 5 wells comb    | 2        |

|                      |   |
|----------------------|---|
| Adjustable slide leg | 4 |
| Electrical Cables    | 1 |